

In order to be confident in your sorting results, cell sorter performance must be evaluated through the appropriate QC methods. **Rmax** is a method that calculates the **maximum recovery** of the sort sample by looking at how much is lost in the unsorted fraction. Deviations from maximum recovery indicate a problem in sorter performance.

Checking post-sort purity is the gold standard to QC sorting performance. However, any issue affecting sorter performance, such as the drop-charge delay timing being off, will affect recovery before affecting purity. Therefore, to assess cell sorter performance, it's important to use a method that measures recovery, rather than purity.

Rmax

During sort, the absolute number of target & non-target particles (t) derived from the **original sample** (Ot/Ont) will end up either in the **sorted** (St/Snt) or the **center stream catch or CSC** (Ct/Cnt), i.e., $O_t = C_t + S_t$ & $O_{nt} = C_{nt} + S_{nt}$

$$R_{max} = \frac{\frac{C_{nt}}{C_t} - \frac{O_{nt}}{O_t}}{\frac{C_{nt}}{C_t} - \frac{S_{nt}}{S_t}}$$


Assuming purity is close to 100% (or above 95%, which it typically is after measuring the drop-charge delay) we can simplify the above equation:

$$R_{max} = 1 - \frac{O_{nt}}{O_t} \times \frac{C_t}{C_{nt}}$$


Rmax can also be used to:

- troubleshoot the instrument
- determine daily drop-charge delay
- evaluate & troubleshoot sample preparation

Rmax online protocol



Calculate Rmax



Riddell (2015) doi: 10.1016/j.jymeth.2015.02.017