




In flow cytometry, the most common data scaling methods are **linear** and **logarithmic**. When presenting data where there is a wide range of fluorescent signal, data is typically displayed in log so all events can be visualized. Often, when analyzing data in log scale, there can be a large portion of events stacked against the axis at zero. Through **biexponential log transformation**, we can visualize the distribution of these negative events around zero in a better way.



Comparison of scaling across same sample

Log: Negative events are against the axis and populations are not well resolved.

Biexponential: Default biex transformation no longer shows events on the axis, but the data is compressed together.

Optimized Biexponential: Clear distribution of all populations, allowing for best visualization and optimal conditions for gating.



Helpful Hints

- Scaling transformation does not alter the data, it simply changes the way the data is visualized.
- Distribution of events around or below zero can be a result of baseline restoration or compensation.
- Biexponential scaling needs to be optimized for better visualization.
- Incorrect scaling of data can hide poorly compensated or unmixed data.
- Scaling adjustments made in acquisition software do not carry over into analysis software.
- Other methods of log transformation can include Hyperlog and ArcSinh.

Useful read: Herzenberg LA, Tung J, Moore WA, Herzenberg LA, Parks DR. Interpreting flow cytometry data: a guide for the perplexed. Nat Immunol. 2006 Jul;7(7):681-5. ([link](#))