



In flow cytometry experiments, **sample preparation is a key component** for having good quality data. The washing steps in a staining protocol allow for the removal of excess antibody. Single washing protocols may not be sufficient to remove this unbound antibody, which can negatively impact the Stain Index.

$$\text{Stain Index} = \frac{\text{median}_{\text{pos}} - \text{median}_{\text{neg}}}{2\sigma_{\text{neg}}}$$




**Fig.1** Stain index graph comparison between **one** and **two** washes in different PE-CF594-CD4 concentrations (N=3).

The Stain Index is higher for samples that have been washed twice compared to those washed only once.



**Fig.2** Graph representation between **one** and **two** washes in different PE-CF594-CD4 concentrations (N=3).



**Fig.3** Graph representation of MFIpos-MFIneg and rSD between **one** and **two** washes in different PE-CF594-CD4 concentrations (N=3).

Experiment above shows in detail the comparison of an antibody titration of CD4 PE-CF594 stained lymphocytes that were washed either one or two times. While the difference between positive and negative Median Fluorescent Intensity (MFI) does not show dramatic changes in the comparison, the Robust Standard Deviation (rSD) of the negative population is significantly lower when the samples were washed twice, resulting in a higher Stain Index. **Including more washing steps in the staining protocol can improve resolution of populations and result in better quality data.**